

DEVELOPMENT OF REAL TIME REACTIVE OXYGEN SPECIES MONITORS

Zoran Ristovski

International Laboratory for Air Quality and Health,
ARC Center for Free Radical Chemistry and Biotechnology,
Biofuels Engine Research Facility



Queensland University of Technology



International Laboratory for Air Quality and Health (ILAQH)



Lidia Morawska

QUT International Laboratory
for Air Quality and Health
Queensland University of Technology

World Health Organization
Collaborating
Centre for
Research and
Training on
Global Burden of
Disease due to
Air Pollution

Contact details:
International Laboratory for Air Quality and Health
Queensland University of Technology
2 George Street, Brisbane Q 4001, Australia
Phone: +61 7 3864 2616/1129
Fax: +61 7 3864 9079
Email: ilaqh@qut.edu.au
Web: <http://www.ilaqh.qut.edu.au>

....since 1993

....until 2003, Environmental
Aerosol Laboratory

....part of IHBI

....since 2004,
WHO Collaborating Centre

....multidisciplinary team of over 25
researchers



a university for the **real** world[®]





Free Radical Chemistry and Biotechnology

An Australian Research Council Centre of Excellence
www.freeradical.org.au



Steven Bottle

Research



Fundamental Radical Chemistry



Good Health & Preventing Disease



Radicals in the Environment



Radicals in Materials Technology



Kinetics & Mechanisms Facility




Pain & Inflammation



Surface Coatings



Climate Change & Energy



Branka Miljevic

5 Universities, >100 Scientists, >\$28 M budget
Specialist MS, Laser, Synthesis, Magnetic Resonance and Fluorescence Equipment



Svetlana Stevanovic



QUT labs: Headed by Prof. Steven Bottle, with Prof. Zoran Ristovski (Associate) and Drs. Miljevic and Fairfull-Smith, 1 PD and 12 PhD students




Biofuels Engine Research Facility - BERF



Richard
Brown

Engine (s)



CONTROL ROOM

CAPABILITIES

- Fuel Handling
 - Liquid and gaseous fuels (Hydrogen)
 - Dual fuel capability
 - Separate Biofuel / Petroleum handling
- Emission Tests
 - Gas (CO₂, CO, HC, NO_x)

Particles




a university for the **real** world[®]



PM & health effects

- Epidemiological studies - strong associations between levels of ambient particulate matter (PM) and increased respiratory and cardiovascular disease morbidity and mortality
- mechanism(s) by which particles induce adverse health effects are still not entirely understood

Proposed mechanism: Oxidative stress hypothesis





PM → free radicals; ROS → oxidative stress → inflammation



cell injury / death

Health effect of PM





Free radicals


- Free radicals can be generated in the human body
- Free radicals make chain reaction
- Antioxidant terminate the chain reaction

IMPLICATED DISEASE STATES




ROS

- Sources: Exogenous and Endogenous
- Oxidative stress



PM & health effects

- **ROS** – Reactive Oxygen Species ($O_2^{\cdot-}$, HO^{\cdot} , RO^{\cdot} , ROO^{\cdot} , 1O_2 , $ONOO^-$, H_2O_2 , $ROOH$)
- **oxidative stress** - an imbalance between the production of ROS and the cell's (or body's) natural antioxidant defence.



Endogenous and exogenous ROS

- ROS can be formed **endogenously**, by the lung tissue cells, during the phagocytic processes initiated by the presence of PM in the lungs, or by particle-related chemical species that have the potential to generate ROS.
- In addition to the particle-induced generation of ROS, several recent studies have shown that particles may also contain ROS (so called, **exogenous ROS**).



Measurements of the radical generation capacity of the PM

Toxicological studies



In vivo


- **The best predictor of PM toxicity**
- *Complicated*
- *Expensive*
- *Time-consuming*
- *Ethical approval*
- *Not suitable for routine analyses and field studies*

In vitro

- **Alternative for in vivo studies**
- *Complicated*
- *Expensive*
- *Time-consuming*
- *Not suitable for routine analyses and field studies*

Acellular assay







- **Cheaper; Faster**
- **No need for ethical approval**
- **Easily applied in field studies and routine analyses**
- *Only inherent oxidative potential; No biological interaction*



DTT




- Dithiotheritol
- Strong reducing agent
- Measure the formation of the ROS by quinones
- Remaining DTT reacts with Ellman reagent
- Following by TNB production in 412 nm
- Report as:
 - nmol DTT min⁻¹μg PM⁻¹
 - normalized index of oxidant generation and toxicity (NIOG)
- *DTT and TNB are both light sensitive*
- *DTT is reactive toward limited number of species*
- *Needs incubation time up to 90 min*



DCFH-DA





- Dichlorofluorescein Diacetate
- Emit fluorescence after being oxidized by hydrogen peroxide
- Detection of H_2O_2
- Needs catalyst usually HRP
- Calibration curve base on hydrogen peroxide equivalent particles
- *Prone to outooxidation*
- *HRP make three fold increase in the fluorescent intensity*




Profluorescent nitroxides

- Suppressed fluorescence emission in the presence of nitroxide moiety
- These molecules react with radicals, leading either to reduction of the nitroxides to the hydroxylamines or oxidation to oxoammonium cation




Profluorescent nitroxides synthesised at QUT

BPEAnit




Phenanthrene-nitroxide

$\lambda_{\text{ex}} = 294 \text{ nm}$
 $\lambda_{\text{em}} = 355 \text{ nm}$
 372 nm




Flourescein-nitroxide

$\lambda_{\text{ex}} = 495 \text{ nm}$
 $\lambda_{\text{em}} = 515 \text{ nm}$



9,10-diphenylanthracene-nitroxide

$\lambda_{\text{ex}} = 395 \text{ nm}$
 $\lambda_{\text{em}} = 410 \text{ nm}$
 430 nm



9,10-bis(phenylethynyl)anthracene-Nitroxide (BPEAnit)

$\lambda_{\text{ex}} = 430 \text{ nm}$
 $\lambda_{\text{em}} = 485 \text{ nm}$
 510 nm

Fairfull-Smith and Bottle. Eur J Org Chem (2008) (32) pp. 5391-5400


Profluorescent nitroxides (PFNs)

9,10-diphenylanthracene-Nitroxide

$\lambda_{ex} = 395 \text{ nm}$

$\lambda_{em} = 410 \text{ nm}$

430 nm



9,10-bis(phenylethynyl)anthracene-Nitroxide (BPEAnit)

$\lambda_{ex} = 430 \text{ nm}$

$\lambda_{em} = 485 \text{ nm}$

510 nm

Fairfull-Smith and Bottle. Eur J Org Chem (2008) (32) pp. 5391-5400

MOTIVATION

Develop a cell-free assay for rapid and routine screenings of the oxidative potential of PM.

Develop a stationary continuous monitor.


Implement the assay in a portable personal exposure sensor.




a university for the **real** world[®]



Particle sampling




BPEAnit assay – sampling:



- bubbling aerosol through an impinger with fritted nozzle tip containing BPEAnit solution (DCFH) fluorescence measurement
- solvent – dimethylsulfoxide (DMSO)
- test & HEPA-filtered control sample taken

- $I_{485\text{nm}}(\text{test}) - I_{485\text{nm}}(\text{ctrl}) \longrightarrow I_{485\text{nm}}(\text{ROS}_{\text{particle}})$
 - ↓ calibration curve
 - ↓ $n(\text{ROS}_{\text{particle}})$
 - ↓ Normalized to the measured particle mass




Real Time Detection of ROS

- Semi-automated systems where particles are collected onto filters first (DTT)
- To skip the extraction procedures particles have to be collected directly into a liquid
 - Water – DCFH
 - DMSO – BPEAnit
- The chemical assay can either be in the collection liquid or can be added later to the collection liquid.




Semi-automated systems


Particle first collected on quartz filters over 23 h. Extraction procedure: 1-inch-diameter punches were extracted in 15 mL of DI water by **sonication**. 5 mL were filtered and used in the DTT assay.




Semi-automated systems



Automated systems based on PILS




DCHF based systems – Clarkson Uni




Venkatachari, P., & Hopke, P. K. (2008)..
Aerosol Science and Technology, 42(8),
629–635

Field results




Wang, Y., Hopke, et al. (2011). Laboratory and Field Testing of an Automated Atmospheric Particle-Bound Reactive Oxygen Species Sampling-Analysis System. *Journal of Toxicology*, 2011(ID 419476)

PFN based system in development - QUT




Orsini, D. A. et al. (2008). A water cyclone to preserve insoluble aerosols in liquid flow - An interface to flow cytometry to detect airborne nucleic acid. *Aerosol Science and Technology*, 42(5), 343–356.

DTT based systems in development - University of Illinois Urbana Champaign




Other non PLS based systems – DCFH (Cambridge Uni)



Fuller, S. J., et al (2014). *Atmospheric Environment*, 92(C), 97–103.

Other non PILS based systems - Impinger based system QUT-HKCityU



Conclusion

- Most systems still in development phase with only 1 of them used as a real time monitor (Clarkson Uni).
- Large problems with auto-oxidation of the chemical probes especially DCFH.
- DTT very complicated to be used.



PERSONAL ROS SENSORS



a university for the **real** world[®]



Approach



Synthesis of novel
Profluorescent
Nitroxide probe
molecules



Cooperative Research Centre for
Polymers
Solutions for a better world

Embedding PFN
probes in polymer
matrices



**International Laboratory
for Air Quality
and Health**
ihbi

Institute of Health and Biomedical Innovation

Implementation of the
sensing technology in a
portable personal sampler



a university for the **real** world[®]




CURRENT TECHNOLOGY

Embedding PFN's into polymers for polymer degradation studies







a university for the **real** world[®]






a university for the **real** world[®]





a university for the **real** world[®]

Imbedding PFN's into polymers



Conclusion

- Technology still in development
- Problems of auto-oxidation of PFN's still not overcome
- Still not sufficient sensitivity for atmospheric levels of ROS.



a university for the **real** world[®]



Acknowledgment

- Organisers for the invitation to present
- Funding:
 - Australian Coal Association Research Program
 - Australian Research Council
 - Institute for Future Environments



a university for the **real** world[®]



Thank you



a university for the **real** world[®]



3
8

Zoran Ristovski's group at ILAQH

