

# PHOTOACOUSTIC spectroscopy utilizing AFM-based CANTILEVER detection

Z. Zelinger, P. Janda

*J. Heyrovský Institute of Physical Chemistry AS CR,  
Dolejškova 3, 182 23 Prague 8 Czech Republic*



The sensitivity in the photoacoustic spectroscopy is limited by the microphone. The replacement of a microphone with a cantilever one enhances remarkably the sensitivity of photoacoustic gas detection:

J. Kauppinen, K. Wilcken, I. Kauppinen, V. Koskinen, "High sensitivity in gas analysis with photoacoustic detection", *Microchem. J.* **76** (2004)151-159.

V. Koskinen, J. Fonsen, J. Kauppinen, I. Kauppinen, "Extremely sensitive trace gas analysis with modern photoacoustic spectroscopy", *Vibr. Spectrosc.* **42** (2006) 239 - 242.


The combination of sensitive microphones and micromechanical elements with laser techniques - new approaches to photoacoustic detection and for chemical sensing possibilities.

graphene leafs:  
outstanding electromechanical properties  
and impressive sensitivity as a mass  
detector  
its utilization as nano/micro-lever sensing  
devices for chemical analysis:


P. Li, Z. You, T. Cui, “Graphene cantilever beams for nano switches“, *Applied Physics Letters*, **101**, (2012) Issue 9, id. 093111.

We have tested:


- AFM-based silicon cantilevers in a square arrangement and
- graphene sheets in the form of a circular membrane.



A044 | 1280x1024 | 2013/11/28 11:13:40




EUROPEAN COOPERATION IN SCIENCE AND TECHNOLOGY



European Network on New Sensing Technologies for Air-Pollution Control and Environmental Sustainability - EuNetAir



# Scheme of detection part consisting silicon cantilevers:



European Network on New Sensing Technologies for Air-Pollution Control and Environmental Sustainability - EuNetAir





Membrane for photoacoustic detector has been prepared from multilayer graphene (MLG) by micromechanical cleavage which is an essential component of sheets stripped-off basal plane highly ordered pyrolytic graphite (HOPG, ZYH Grade, Bruker, USA) :

Novoselov, K. S.; Geim, A. K.; Morozov, S. V.; Jiang, D.; Zhang, Y.; Dubonos, S. V.; Grigorieva, I. V.; Firsov, A. A. (2004). "Electric Field Effect in Atomically Thin Carbon Films". *Science* 306 (5696): 666–669.



MLG leafs (thickness  $<10 \mu\text{m}$ ) have been mounted on glass window of cuvette for photoacoustic spectroscopy and used for detection of acoustic waves by laser beam reflected to CCD detector.




European Network on New Sensing Technologies for Air-Pollution Control and Environmental Sustainability - EuNetAir

The logo for EuNetAir features the text "EuNetAir" in a stylized font, with "u" and "et" in smaller letters. To the left is the European Union flag, and to the right is a mobile phone with a signal icon. Below the logo is a globe and a chemical structure. The text "EUUNETAIR" is written in a bold, sans-serif font below the globe and chemical structure.

European Network on New Sensing Technologies for Air-Pollution Control and Environmental Sustainability


# Scheme of detection part based on graphene:




European Network on New Sensing Technologies for Air-Pollution Control and Environmental Sustainability - EuNetAir

European Network on New Sensing Technologies for Air-Pollution Control and Environmental Sustainability


# Experimental set-up



# Detected MLG leafs movements




4 Hz





23 Hz

# Comparison of the same concentration detection by microphone and graphene processed by lock-ins – laser radiation is ON or OFF



Sensitivity levels of the investigated experimental set-ups have been tested by the utilizing of the concentration standards based on the **permeation method** in the flow regime.

# ON and OFF concentration (~2 ppm) for microphone and graphene detection:





- silicon cantilevers and multilayer graphene (MLG) leafs as photoacoustic detectors.
- future:  
few-layer graphene (FLG)  
cantilevers/membranes, prepared by CVD  
technique and transferred to support,  
allowing free-standing mounting.

Many thanks for your attention!